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Email : bip@clickTA.com

PLASMA OR SERUM ?

Serum is obtained from whole blood by centrifugation after completion of the platelet and clotting factor coagulation process. Serum must therefore be regarded as an artifact. By definition, it is devoid of clotting factors but is enriched with the cellular components of platelets and metabolic products.

Plasma is the virtually cell-free supernatant following centrifugation of whole blood, the coagulability of which is inhibited by the addition of anticoagulants immediately after sampling. Anticoagulants inhibit clot formation through various mechanisms.

For some analytes, there are diagnostically relevant differences between the results obtained from serum and those obtained from plasma (see Table 1 and Figure 1)


Table 1 Analytes with diagnostically relevant serum / heparinized plasma concentration differences and their main causes (1)

Analyte	% change in comparison to the mean in plasma	Main cause of the serum / plasma difference
Potassium	+6.2	Lysis of the cells, particularly the platelets*
Inorganic phosphate	+10.7	Release from cellular elements
Total protein	-5.2	Effect of fibrinogen
Ammonia	+38	Thrombocytolysis, hydrolysis of glutamine
Lactate	+22	Release from cellular elements

*Also causally responsible for high potassium values in serum are – in addition to thrombocytolysis – haemolysis and extreme leucocytolysis (when leucocyte counts and > 50 G/L). Concerning the platelet influence in whole blood, the following applies an increase in the platelet count by 100 G/L means an average increase of 0.11 mmol/L in the difference between the serum and plasma potassium values.


Figure 1 Plasma serum differences obtained in 4 ml separation tubes (1). Ratio of the median difference between serum and plasma and the coefficient of variation of the analytical procedure used.

A linear correlation has been shown between the difference in serum-plasma potassium and phosphate and the number of platelets in blood (Fig 2). Some constituents found in high concentrations in platelets can therefore not be correctly determined in serum, e.g. acid phosphatase activity, neuron specific enolase, dopamine and serotonin (2).

Differences between the values obtained in serum and plasma are due to the following physiological and technical reasons:

- The analyte may be used up during clotting: fibrinogen, platelets, glucose.
- The analyte may be released from cells during clotting : potassium, lactate dehydrogenase, phosphate, ammonia, lactate.
- The anticoagulant may interfere with the assay or contaminate some assays: Y-glutamyl transferase, lithium in flame photometry, when calibrated with lithium.
- Methodology of the respective determination, including for example the type of measuring instrument used (monochromatic/bichromatic measurement) (4) may cause interference.
- Fibrinogen may interfere with the method (some heterogeneous immuno assays).


Figure 2

Dependence of plasma –serum difference in potassium on platelet count in blood (3).

What advantages does plasma have over serum?

1. **Time saving.** Waiting for blood to clot is eliminated. The centrifugation period can be reduced considerably by increasing the rotation speed.
2. **Higher yield.** Approx. 15-20 % more plasma than serum can be obtained from whole blood.
3. Virtually no interference due to subsequent coagulation. Post-centrifugal coagulation can occur in serum. This effect does not occur in plasma.
4. Results from plasma are more representative of the in vivo state compared to serum.
4. Low risk of hemolysis and thrombocytolysis. In healthy persons, free hemoglobin is about 10 times less concentrated in plasma than in serum. In plasma, the platelets remain intact in vitro; there is hence no pseudohyperkalemia here, as is found in serum.

What are the disadvantages of plasma relative to serum?

1. Protein electrophoresis is altered. Fibrinogen appears as a peak in the region of the Y-globulins and can simulate or mask an M-gradient.
2. Method-dependent interference. Anticoagulants can - as potential complexing agents and enzyme inhibitors-lead to method-dependent interference. Every new procedure should therefore be tested for , eg., heparin interference (**see leaflets**)
3. Cation interference. When heparinates are used, lithium or ammonium may interfere with the methods for determining them.

Blood Collection Reminders

- ◆ Disposable blood sampling systems with additives facilitate and improve the sampling procedure and have made a world – wide contribution to a certain degree of standardization in the collection of venous blood.
- ◆ If serum of plasma separator tubes are used, renewed centrifugation following storage of the sample in the refrigerator should be avoided as this would lead to noticeable increases in the concentration of potassium, inorganic phosphate, lactate dehydrogenase and alanine aminotransferase activities.
- ◆ When obtaining a blood sample, it is imperative that attention be paid to ensuring thorough mixing of the blood with the anticoagulant; foaming should be avoided.
- ◆ Not more than 2 minutes should elapse between the beginning of stasis and the mixing of the blood with anticoagulant.

- ◆ When obtaining serum, it is necessary to wait until the blood has clotted completely (30 min at room temperature); if this is not done, post-coagulation, which can cause volume errors in dispensing, **as well as blockages in the analyzer system can occur.**

Types of plasma : For coagulation testing various types of plasma obtained from citrated blood are used:

Table 2 Different types of plasma

Plasma	Relative centrifugal force (g)	Centrifugation time (min)
Platelet-rich	150-200	5
Platelet-poor	1,000-2,000	10
Platelet-free	2,000-3,000	15-30

References

1. Voit R. Plasma-Serum-Unterschiede und lagerungs stabilität. Kinisch-chemische Meßgrößen bei verwendung von Plasmatrennröhrchen. Dissertation, Ludwig-Maximilians-Universität Munchenc, 1993.
2. Guder Wg. Thrombocytes as interfering factors in Clinical Chemistry. Editorial Comment. J Clin Chem Clin Biochem 1990; 28: 445.
3. Lutomski DM. Bower RH. The effect of thrombo cytosis on serum potassium and phosphorus concentration. Am J Med Sci 1994; 307:255-8.
4. Hallbach J, Hoffmann GE, Guder WG. Overestimation of albumin in heparinized plasma. Clin Chem 1991; 37:556-86.
5. W.G. Guder, S. Nrayanan, Smaple From the Patient to the Laboratory, 1996

PHLEBOTOMIST AND LABORATORY'S UNIVERSAL PRECAUTIONS.

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